


Product Data Sheet

Sight Glass LWZ 140, LWZ 141

Accessory for Integral



Manufacturer:

LAUDA DR. R. WOBSE GMBH & CO. KG

Laudaplatz 1

97922 Lauda-Königshofen

Germany

Telephone: +49 (0)9343 503-0

E-mail: info@lauda.de

Internet: <https://www.lauda.de>

Translation of the original

Q4DA-E_17-007, 1, en_US 5/22/2024 © LAUDA 2024

Product Data Sheet Sight Glass Stainless Steel




Fig. 1: Sight Glass

Technical Data

Sight glass stainless steel LWZ 140
Connection: M30 x 1,5

Sight glass stainless steel LWZ 141
Connection: M38 x 1,5

Data	Value	Unit
Permissible pressure range PS	-0,5 – 12	bar
Operating temperature range TS	-40 – 140	°C



The sight glass must not be operated outside the approved ranges.

Manufacturer:

LAUDA DR. R. WOBSE GMBH & CO. KG • Laudaplatz 1 • 97922 Lauda-Königshofen

Telephone: +49 (0)9343 503-0

E-mail: info@lauda.de • Internet: <https://www.lauda.de>