

## PRODUCT DATA SHEET



### VS 15 T

LAUDA Varioshake shakers execute a tumbling movement of 2 to 50 rpm for reproducible preparation and ideal growth conditions e.g. for cell cultures in bags, petri dishes or edge plates.

Technical Attributes	VS 15 T
Ambient temperature min.	10 °C
Ambient temperature max.	50 °C
Load capacity max.	15 kg
Shaking amplitude	3 Winkelgrade zur Waagerechten
Shaking motion	three-dimensional, orbital-rocking
Timer	1 minute - 99:59 hours or continuous operation
Dimensions (WxDxH) in mm	480x487x160
Weight	19.5 kg

LAUDA DR. R. WOBSEY GMBH & CO. KG  
Pfarrstraße 41/43 · 97922 Lauda-Königshofen  
Postfach 1251 · 97912 Lauda-Königshofen · DE

T +49 (0) 9343 503-0 · F +49 (0) 9343 503-222  
info@lauda.de · www.lauda.de  
WEEE-Reg.-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen  
Registergericht Mannheim · HRA 560069

Persönlich haftende Gesellschafterin:  
LAUDA DR. R. WOBSEY Verwaltungs-GmbH  
Sitz Lauda-Königshofen  
Registergericht Mannheim · HRB 560226

Geschäftsführer:  
Dr. Gunther Wobser (Vors.), Dr. Mario Englert  
Dr. Marc Stricker  
Beirat: Dr. Gerhard Wobser