

## PRODUCT DATA SHEET



### VS 15 R

The robust low-noise rocking shakers from the Varioshake product line are sure to impress with a reduced footprint in the laboratory as well as load-independent speed and digital control. They are perfect for continuous operation.

Technical Attributes	VS 15 R
Ambient temperature min.	10 °C
Ambient temperature max.	50 °C
Shaking platform width	450 mm
Shaking platform depth	450 mm
Load capacity max.	15 kg
Shaking amplitude	3 Winkelgrade zur Waagrechten
Shaking frequency min.	2 min <sup>-1</sup>
Shaking frequency max.	50 min <sup>-1</sup>
Shaking motion	rocking
Timer	1 minute - 99:59 hours or continuous operation
Dimensions (WxDxH) in mm	480x487x160
Weight	19.5 kg

LAUDA DR. R. WOBSEY GMBH & CO. KG  
Pfarrstraße 41/43 · 97922 Lauda-Königshofen  
Postfach 1251 · 97912 Lauda-Königshofen · DE

T +49 (0) 9343 503-0 · F +49 (0) 9343 503-222  
info@lauda.de · www.lauda.de  
WEEE-Reg-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen  
Registergericht Mannheim · HRA 560069

Persönlich haftende Gesellschafterin:  
LAUDA DR. R. WOBSEY Verwaltungs-GmbH  
Sitz Lauda-Königshofen  
Registergericht Mannheim · HRB 560226

Geschäftsführer:  
Dr. Gunther Wobser (Vors.), Dr. Mario Englert  
Dr. Marc Stricker  
Beirat: Dr. Gerhard Wobser