


PRODUCT DATA SHEET

Versafreeze VF 75085 C

LAUDA Versafreeze chest freezers offer state-of-the-art refrigeration technology with a temperature range of 0 to -85 °C, excellent insulation properties and the greatest possible protection thanks to password-protected access rights and secure monitoring via the LAUDA Cloud. LAUDA Versafreeze freezers are the safe storage place for your valuable samples and pharmaceuticals.



Working temperature min.
-85 °C



Working temperature max.
-50 °C

Technical Attributes

Working temperature min.	-85 °C
Working temperature max.	-50 °C
Ambient temperature min.	15 °C
Ambient temperature max.	28 °C
Usable volume	754 L
Interior dimensions (WxDxH) in mm	1600 x 620 x 760
Cool down time	from 20 °C to -85 °C in 5 h
Thawing time	from -85 °C to 0 °C in 19 h
Dimensions (WxDxH) in mm	2102 x 910 x 1056

Versafreeze VF 75085 C

LAUDA DR. R. WOBSEY GMBH & CO. KG
Pfarrstraße 41/43 • 97922 Lauda-Königshofen
Postfach 1251 • 97912 Lauda-Königshofen • DE

T +49 (0) 9343 503-0 • F +49 (0) 9343 503-222
info@lauda.de • www.lauda.de
WEEE-Reg.-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert
Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser