

## PRODUCT DATA SHEET

### Proline PB C

LAUDA Proline bridge thermostats are available in two versions with different pump models and immersion depths. The PB models have a pressure/suction pump and require a bath depth of 200 mm, while the PBD models have a more powerful pressure pump (D) and require a bath with a depth of 320 mm. In addition, both series of models differ in the selected control head: Master or Command (C). Through variably extendable telescopic rods, all models can be attached without problem to baths with a width from 310 mm up to 550 mm.



Working temperature min.  
30 °C



Working temperature max.  
300 °C

#### Technical Attributes

Working temperature min.

Working temperature max.

Ambient temperature min.

Ambient temperature max.

Temperature stability

Application

Heater power

Power consumption

Dimensions (WxDxH) in mm

Pump pressure max.

Pump flow max. (pressure)

Pump connection thread

#### Proline PB C

30 °C

300 °C

5 °C

40 °C

0.01 ±K

external | Inside the bath | External and inside the bath

3.5 kW

3.6 kW

-x185x520

0.7 bar

25 L/min

M16 x 1

LAUDA DR. R. WOBSEY GMBH & CO. KG  
Pfarrstraße 41/43 · 97922 Lauda-Königshofen  
Postfach 1251 · 97912 Lauda-Königshofen · DE

T +49 (0) 9343 503-0 · F +49 (0) 9343 503-222  
info@lauda.de · www.lauda.de  
WEEE-Reg.-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen  
Registergericht Mannheim · HRA 560069

Persönlich haftende Gesellschafterin:  
LAUDA DR. R. WOBSEY Verwaltungs-GmbH  
Sitz Lauda-Königshofen  
Registergericht Mannheim · HRB 560226

Geschäftsführer:  
Dr. Gunther Wobser (Vors.), Dr. Mario Englert  
Dr. Marc Stricker  
Beirat: Dr. Gerhard Wobser