

PRODUCT DATA SHEET

Standing: 2023-07-06

Pt100 RTD (resistance temperatur detector)

Measuring range: -100 °C ... 200 °C

Part Number: ETP 059

Features

- Push-in resistance thermometer with permanently connected silicone cable (length: 2 m), silicone cable shielded
- Connector plug mounted directly on the probe.
- Connector plug: Lemo connector, size 1, four-pole
- Resistance thermometer (4-wire circuit) class A according to DIN EC 751



Technical Features

Usable floor space width	300 mm
Material	CrNi steel (1.4571)
Outer diameter	6 mm
Immersion depth	250 mm

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0 • F + 49 (0) 9343 503-222
info@lauda.de • www.lauda.de
WEEE-Reg.-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Ralf Hermann, Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser