

PRODUCT DATA SHEET

Standing: 2023-07-09

Stainless steel inset rack
for 72 eppendorf tubes
Part Number: UE 028

Features

- Placing the rack on the edge of the bath.
- Hole grid: 5 mm, hole diameter: 19 mm
- Perforation: 7 x 6 holes



Technical Features

Overall dimensions (WxDxH)	154 x 100 x 21.5 mm
Usable floor space width	134 mm
Usable floor space depth	80 mm
Material	Stainless steel
Outer diameter	11 mm
Number of samples	42
Immersion depth	21.5 mm

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0 • F + 49 (0) 9343 503-222
info@lauda.de • www.lauda.de
WEEE-Reg.-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Ralf Hermann, Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser