

PRODUCT DATA SHEET

Standing: 2023-07-09

Stainless steel inset rack RN 30
for 15 tubes, 24-30 mm Dia. 110 mm ID
Part Number: UE 032

Features

- Rack with plastic handle(s) for a temperature up to 100 °C
- Distance between levels: 50 mm
- Perforation: 3 x 5 holes



Technical Features

Overall dimensions (WxDxH)	300 x 145 x 137 mm
Material	Stainless steel
Number of samples	15
Immersion depth	110 mm

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0 • F + 49 (0) 9343 503-222
info@lauda.de • www.lauda.de
WEEE-Reg.-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Ralf Hermann, Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser