

PRODUCT DATA SHEET

Standing: 2023-07-09

Stainless steel inset rack RQ 18/1
for 21 tubes, Ø 14-18 mm 80 mm ID
Part Number: UG 106

Features

- For comfortable handling in bath thermostats and water baths.
- Rack with plastic handle(s) for a temperature up to 100 °C
- suitable for test tubes with a diameter of 14 - 18 mm
- Distance from bottom: 62 mm

Technical Features

Overall dimensions (WxDxH)	190 x 95 x 183 mm
Material	Stainless steel
Number of samples	21
Immersion depth	80 mm

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0 • F + 49 (0) 9343 503-222
info@lauda.de • www.lauda.de
WEEE-Reg.-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Ralf Hermann, Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser