

## PRODUCT DATA SHEET

Standing: 2023-07-09

Stainless steel inset rack RQ 30  
for 10 tubes, Ø 24-30 mm 110 mm ID  
Part Number: UG 108

### Features

- For comfortable handling in bath thermostats and water baths.
- Rack with plastic handle(s) for a temperature up to 100 °C
- suitable for test tubes with a diameter of 24 - 30 mm
- Distance from bottom: 60 mm

### Technical Features

Overall dimensions (WxDxH)

190 x 95 x 203 mm

Material

Stainless steel

Number of samples

10

Immersion depth

110 mm

LAUDA DR. R. WOBSEY GMBH & CO. KG  
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0 • F + 49 (0) 9343 503-222  
info@lauda.de • www.lauda.de  
WEEE-Reg.-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen  
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:  
LAUDA DR. R. WOBSEY Verwaltungs-GmbH  
Sitz Lauda-Königshofen  
Registergericht Mannheim • HRB 560226

Geschäftsführer:  
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,  
Dr. Ralf Hermann, Dr. Marc Stricker  
Beirat: Dr. Gerhard Wobser