


PRODUCT DATA SHEET

Published: 2026-08-19

Insert rack
with 372 openings, Ø 13 mm
Part Number: A000020

Features

- Stainless steel test tube rack
- Easy positioning and fixing on the four stainless steel pins in the shaker basket.
- Two handles extending from the water surface allow for easy placement and removal.
- Height of the handle: 100 mm
- Distance from bottom: 54 mm



Technical Features

Overall dimensions (WxDxH)

450 x 300 x 155 mm

Material

Stainless steel

Weight

1.60 kg

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
info@lauda.de • www.lauda.de
WEEE-Reg-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser