

PRODUCT DATA SHEET

Standing: 2024-03-22

milk bottle rack with 12 openings Ø56 mm

Part Number: A000022

Features

- Stainless steel test tube rack
- Two handles extending from the water surface allow for easy placement and removal.
- Height of the handle: 108 mm
- Distance from bottom: 54 mm
- Perforation: 5 x 4 holes



Technical Features

Overall dimensions (WxDxH)	240 x 195 x 185 mm
Material	Stainless steel
Outer diameter	56 mm
Number of samples	12

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
info@lauda.de • www.lauda.de
WEEE-Reg.-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Ralf Hermann, Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser