

PRODUCT DATA SHEET

Published: 2026-06-29

Shaking tray
for Shaking Water Bath
Part Number: A000023

Features

- Shaking tray made of stainless steel.
- With perforated grid to accommodate clamps for Erlenmeyer flasks and to fix test tube racks A000032 and A000033.
- Perforation: 33 x 21 holes
- Hole grid: 12 mm, hole diameter: 4.2 mm
- Two handles extending from the water surface allow for easy placement and removal.
- When using the tray, the useable depth is reduced by 15 mm: Usable depth: 175 mm



Technical Features

Overall dimensions (WxDxH)	450 x 300 x 10 mm
Usable floor space width	400 mm
Usable floor space depth	256 mm
Material	Stainless steel
Weight	0.99 kg

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
info@lauda.de • www.lauda.de
WEEE-Reg.-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser