

PRODUCT DATA SHEET

Published: 2025-09-16

Adhesive mat for shaking platform or tray

200 x 200 mm, can be cut to size

Part Number: A000041

Features

- A universal solution to fasten most various vessels onto the shaking platform or tray
- Max. shaking speed up to 250 rpm, depending on shape and weight of the vessel as well as on the shaking amplitude
- Operation temperature range 15 to 60 °C

Technical Features

Overall dimensions (WxDxH)

200 x 200 x null mm

Ambient temperature range

15 ... 60 °C

Weight

0.21 kg

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
info@lauda.de • www.lauda.de
WEEE-Reg-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser