

## PRODUCT DATA SHEET

Standing: 2024-03-22

Shaking tray for Shaker VS 30 ○  
and Shaking incubator VS 150 ○I  
Part Number: A000047

### Features

- Shaking tray made of anodized aluminium.
- Tray is equipped with four stand rods (length: 150 mm) for secure attachment of additional equipment.
- With perforated grid to accommodate clamps for Erlenmeyer flasks and to fix test tube racks A000059 and A000060.
- Perforation: 11 x 9 holes
- Hole grid: 55 mm, hole diameter: 4.2 mm



### Technical Features

Overall dimensions (WxDxH)	670 x 537 x null mm
Usable floor space width	554 mm
Usable floor space depth	447 mm

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG  
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0  
info@lauda.de • www.lauda.de  
WEEE-Reg.-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen  
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:  
LAUDA DR. R. WOBSEY Verwaltungs-GmbH  
Sitz Lauda-Königshofen  
Registergericht Mannheim • HRB 560226

Geschäftsführer:  
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,  
Dr. Ralf Hermann, Dr. Marc Stricker  
Beirat: Dr. Gerhard Wobser