

PRODUCT DATA SHEET

Published: 2025-12-03

Shaking tray
for Shaker and Shaking Incubator
Part Number: A000047

Features

- Shaking tray made of anodized aluminium.
- Tray is equipped with four stand rods (length: 150 mm) for secure attachment of additional equipment.
- With perforated grid to accommodate clamps for Erlenmeyer flasks and to fix test tube racks A000059 and A000060.
- Perforation: 11 x 9 holes
- Hole grid: 55 mm, hole diameter: 4.2 mm



Technical Features

Overall dimensions (WxDxH)

670 x 537 x 3 mm

Usable floor space width

554 mm

Usable floor space depth

447 mm

Weight

2.90 kg

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
info@lauda.de • www.lauda.de
WEEE-Reg-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser