

PRODUCT DATA SHEET

Published: 2026-06-29


Insert rack for 24 test tubes

Length: 75 - 160 mm, Ø 12-17 mm

Part Number: A000059

Features

- suitable for test tubes Ø12 - 17 mm, e.g. 15 mL falcon tubes
- The perforated part can be swiveled by 90 ° and is equipped with pressure springs for a tight hold.
- The rack must be screwed onto a tray, e.g., one of the Varioshake series shaking trays A000044 to A000047. Not suitable for use in shaking water baths.
- Stainless steel test tube rack



Technical Features

Overall dimensions (WxDxH)

340 x 60 x 80 mm

Number of samples

24

Weight

0.61 kg

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
info@lauda.de • www.lauda.de
WEEE-Reg.-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser