

PRODUCT DATA SHEET

Published: 2025-09-16

Rack for 16 tubes

Ø 25-29 mm

Part Number: A000060

Features

- suitable for test tubes Ø25 - 29 mm, e.g. 50 mL falcon tubes
- The perforated part can be swiveled by 90 ° and is equipped with pressure springs for a tight hold.
- The racks must be screwed onto a tray, e.g. shaking trays A000044 to A000047 for Varioshake. Not suitable for the use in shaking water baths.
- Stainless steel test tube rack



Technical Features

Overall dimensions (WxDxH)

340 x 60 x 80 mm

Number of samples

16

Weight

0.77 kg

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
info@lauda.de • www.lauda.de
WEEE-Reg-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser