

PRODUCT DATA SHEET

Published: 2026-06-29

Temperature brake
for data logger LOG100CRYO (A001383)
Part Number: A001384

Features

- Temperature retarder made of aluminum for installation in freezers
- A monitoring sensor can be attached directly in one of the holes in the temperature retarder. This prevents premature alarms when the door/lid is opened briefly
- Three holes Ø 3.5mm / 4.2mm / 5mm, each with a hole for a locking screw



Technical Features

Overall dimensions (WxDxH)

60 x 25 x 25 mm

Weight

0.20 kg

LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
info@lauda.de • www.lauda.de
WEEE-Reg-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser