

PRODUCT DATA SHEET

Published: 2025-09-16

Storage rack for chest freezers 556 l und 754 l
for 37 test plates
Part Number: A001400

Features

- Thermostatic bath/circulator with simplest operation and integrated cooling system
- Proprietary CAN protocol based on two freely configurable messages with multiplexed signals
- simple and uncomplicated for temporary use in projects, for demo and testing processes

Technical Features

Overall dimensions (WxDxH)	91 x 140 x 746 mm
Material	Stainless steel
Weight	1.40 kg

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
info@lauda.de • www.lauda.de
WEEE-Reg-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser