

PRODUCT DATA SHEET

Published: 2025-09-16

Bath cover (Universa)

for bath opening 200 x 200 mm, with 5 ring inserts

Part Number: A001745

Features

- With the bath lid, tall beakers with a nominal volume of between 50 mL and 250 mL can be suspended.
- Easy positioning of the bath lid thanks to the two bow handles. No tools required!
- The individual ring sets are in four parts and made of heat-resistant plastic
- Individual ring sets are removable
- Two additional sockets for external Pt100 temperature sensors (D = 4 mm and 1.9 mm)
- suitable for U 12, U 20, U 1225, U 1245, U 2040

Technical Features

Overall dimensions (WxDxH)

232 x 232 x 47 mm

Material

Stainless steel

Ring Nr.	Durchmesser der Öffnung	passend für ein Becherglas	
		mit einem Nennvolumen	mit einem Außendurchmesser
ohne Ring	61 mm	250 mL	60 mm
1	55 mm	150 mL	54 mm
2	49 mm	100 mL	48 mm
3	39 mm	50 mL	38 mm
4*	16 mm	Reagenzglas, Durchmesser 16 mm	

*Der Gewindeknauf ist abschraubbar, so dass ein Reagenzglas (D = 16 mm) eingesetzt werden kann.

LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
info@lauda.de • www.lauda.de
WEEE-Reg.-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser