

PRODUCT DATA SHEET

Published: 2025-12-03

Stainless steel insert rack
for 30 samples, D = 17 mm
Part Number: A001788

Features

- Ideally suited for Falcon Tubes 15 mL
- Excellent clamping action ensures secure fixation.
- Due to the hanging position, the tubes are evenly rinsed and tempered.
- The retaining bracket is used to load and remove the samples.
- Suitable for bathroom openings measuring 150x130 mm



Technical Features

Overall dimensions (WxDxH)

186 x 128 x 128 mm

Number of samples

30

Material

Stainless steel



LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
info@lauda.de • www.lauda.de
WEEE-Reg-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser