

## PRODUCT DATA SHEET

Standing: 2025-03-18

Pt100 RTD (resistance temperatur detector)

Measuring range: -100 °C ... 200 °C

Part Number: ETP 059

### Features

- Push-in resistance thermometer with permanently connected silicone cable (length: 2 m), silicone cable shielded
- Cable length: 2 m
- Connector plug: Lemo connector, size 1, four-pole
- Resistance thermometer (4-wire circuit) class A according to DIN EC 751



### Technical Features

|                             |                     |
|-----------------------------|---------------------|
| Outer diameter              | 4 mm                |
| Immersion depth             | 210 mm              |
| Operating temperature range | -100 ... 200 °C     |
| Material                    | CrNi steel (1.4571) |
| Length                      | 25 cm               |
| Weight                      | 0.10 kg             |

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG  
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0  
info@lauda.de • www.lauda.de  
WEEE-Reg-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen  
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:  
LAUDA DR. R. WOBSEY Verwaltungs-GmbH  
Sitz Lauda-Königshofen  
Registergericht Mannheim • HRB 560226

Geschäftsführer:  
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,  
Dr. Marc Stricker  
Beirat: Dr. Gerhard Wobser