

PRODUCT DATA SHEET

Published: 2025-09-16

Pt 100 Platinum resistance thermometer
measuring range: -200 °C ... 450 °C, cable: 1.5 m
Part Number: ETP 065

Features

- Immersion probe, sensor class 1/3 DIN
- With black handle and cable 2 m



Technical Features

Outer diameter	3 mm
Immersion depth	300 mm
Material	CrNi steel (1.4571)
Weight	0.18 kg

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
info@lauda.de • www.lauda.de
WEEE-Reg-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser