

PRODUCT DATA SHEET

Standing: 2025-08-04

Platinum resistance therm. Pt100-91

Measuring range: -100 °C ... 300 °C

Part Number: ETP 127

Features

- Connector plug mounted directly on the probe.
- Probe design: straight
- Connector plug: Lemo connector, size 1, four-pole
- Cable resistance thermometer (4-wire circuit) class AA according to EN 60751

Technical Features

Outer diameter	5 mm
Immersion depth	80 mm
Operating temperature range	-100 ... 300 °C
Material	CrNi steel (1.4571)
Weight	0.03 kg

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
info@lauda.de • www.lauda.de
WEEE-Reg-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser