

PRODUCT DATA SHEET

Published: 2026-06-29

Automatic refill
with pump, 24 V, 2x hose nozzles 1/2"
Part Number: LWZ 142

Features

- Refill system for FD 50 with pump for automatic compensation of fluid losses when changing consumers
- Permissible temperature control media: Water/glycol mixtures
- Filling from storage tank by direct suction
- With bracket for the frame profile of the FD 50



Technical Features

Overall dimensions (WxDxH)	150 x 145 x 225 mm
Operating temperature range	10 ... 40 °C
Pump flow rate max. (pressure)	7 L/min
Connection 1 thread (outer)	1/2"
Connection 2 thread (outside)	1/2"

LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
info@lauda.de • www.lauda.de
WEEE-Reg.-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser