

## PRODUCT DATA SHEET

Standing: 2025-01-16

Metal tube MTK200 G3/4" to G3/4"  
DN 20, length 2000 mm  
Part Number: LZM 076

### Features

- Flat sealing



### Technical Features

Operating temperature range

-50 ... 150 °C

Connection 1 thread (inside)

G 3/4"

Connection 2 thread (inside)

G 3/4"

Length

200 cm

Material

Stainless steel

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG  
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0  
info@lauda.de • www.lauda.de  
WEEE-Reg.-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen  
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:  
LAUDA DR. R. WOBSEY Verwaltungs-GmbH  
Sitz Lauda-Königshofen  
Registergericht Mannheim • HRB 560226

Geschäftsführer:  
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,  
Dr. Marc Stricker  
Beirat: Dr. Gerhard Wobser