

PRODUCT DATA SHEET

Published: 2026-07-23

Viton tube 11 mm i.D., insulated

Part Number: RKJ 091

Features

- Viton machine hose, smooth on the inside and outside, density approx. 1.87 g/cm^3 , tensile strength $>7 \text{ N/mm}^2$, elongation $>260 \%$
- Excellent resistance to ozone, weathering, light and ageing: Resistant to oils, most mineral acids and many aliphatic and aromatic hydrocarbons
- The elastomer complies with EC Directive 2002/95/EC and/or 2011/65/EU (RoHS)
- Free of PFOS, free of PTFE and silicone oils
- Hose available by the metre



Technical Features

Operating temperature range

-18 ... 230 °C

Inner diameters

11 mm

Material

FKM

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
info@lauda.de • www.lauda.de
WEEE-Reg-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser