

PRODUCT DATA SHEET

Standing: 2024-11-15

Polycarbonate inset rack RK 075
for 30 tubes, 13/14 mm dia.
Part Number: UE 023

Features

- Use up to $T_{max} = 100\text{ °C}$
- Placing the rack on the edge of the bath.
- Perforation: 6 x 5 holes



Technical Features

Overall dimensions (WxDxH)	184 x 100 x 92 mm
Usable floor space width	108 mm
Usable floor space depth	100 mm
Immersion depth	63 mm
Material	Polycarbonate
Weight	0.2 kg

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
info@lauda.de • www.lauda.de
WEEE-Reg-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser