

PRODUCT DATA SHEET

Published: 2026-06-29

Insert rack, polycarbonate
for 30 tubes, Ø 13-14 mm
Part Number: UE 023

Features

- Use up to $T_{max} = 100\text{ °C}$
- Placing the rack on the edge of the bath.
- Perforation: 6 x 5 holes



Technical Features

Overall dimensions (WxDxH)

184 x 100 x 92 mm

Usable floor space width

108 mm

Usable floor space depth

100 mm

Immersion depth

63 mm

Material

Polycarbonate

Weight

0.16 kg

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
info@lauda.de • www.lauda.de
WEEE-Reg-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser