

PRODUCT DATA SHEET

Standing: 2024-11-15

Rack for 20 tubes

Part Number: UE 049

Features

- Suitable for centrifuge tubes with a diameter of 42 mm to 48 mm
- For use at T = 0 °C ... 135 °C
- Perforation: 4 x 5 holes



Technical Features

Overall dimensions (WxDxH)

250 x 100 x 65 mm

Material

Stainless steel

Weight

1.5 kg

LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
info@lauda.de • www.lauda.de
WEEE-Reg-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser