

PRODUCT DATA SHEET

Published: 2026-06-29

Insert rack
for 10 tubes, Ø 24-30 mm
Part Number: UG 108

Features

- For comfortable handling in bath thermostats and water baths.
- Rack with plastic handle(s) for a temperature up to 100 °C
- suitable for test tubes with a diameter of 24 - 30 mm
- Distance from bottom: 60 mm

Technical Features

Overall dimensions (WxDxH)

190 x 95 x 203 mm

Number of samples

10

Immersion depth

110 mm

Material

Stainless steel

Weight

0.36 kg

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
info@lauda.de • www.lauda.de
WEEE-Reg-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser