

PRODUCT DATA SHEET

Published: 2025-09-16

Insert rack (Rotorack)
with 18 holes, Ø 11 mm
Part Number: UG 110

Features

- suitable for REJ 1225 G
- Accommodates thermometers or temperature probes for calibration purposes.
- Rotorack for 18 specimens with a diameter of 11 mm.
- The insert can be rotated and adjusted in height by 40 mm.



Technical Features

Outer diameter	175 mm
Number of samples	18
Material	Stainless steel
Weight	1.00 kg

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
info@lauda.de • www.lauda.de
WEEE-Reg-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser