

## PRODUCT DATA SHEET

Published: 2026-06-29

Insert rack (Rotorack)  
with 18 holes, Ø 11 mm  
Part Number: UG 110

### Features

- suitable for REJ 1225 G
- Accommodates thermometers or temperature probes for calibration purposes.
- Rotorack for 18 specimens with a diameter of 11 mm.
- The insert can be rotated and adjusted in height by 40 mm.



### Technical Features

Outer diameter	175 mm
Number of samples	18
Material	Stainless steel
Weight	1.00 kg

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG  
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0  
info@lauda.de • www.lauda.de  
WEEE-Reg-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen  
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:  
LAUDA DR. R. WOBSEY Verwaltungs-GmbH  
Sitz Lauda-Königshofen  
Registergericht Mannheim • HRB 560226

Geschäftsführer:  
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,  
Dr. Marc Stricker  
Beirat: Dr. Gerhard Wobser