

PRODUCT DATA SHEET

Published: 2026-06-29

Insert rack (Rotorack)
with 15 holes, Ø 12 mm
Part Number: UG 112

Features

- Accommodates thermometers or temperature probes for calibration purposes.
- Rotorack for 15 specimens with a diameter of 12 mm.
- The insert can be rotated and adjusted in height by 60 mm.



Technical Features

Outer diameter	146 mm
Number of samples	15
Material	Stainless steel
Weight	1.00 kg

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
info@lauda.de • www.lauda.de
WEEE-Reg.-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser